BioRoot™ Flow

Active Biosilicate Technology™

| | |
|---|---------|
| fr - Ciment de scellement canalaire minéral bioactif..... | 2 - 3 |
| en - Bioactive mineral root canal sealer | 4 - 5 |
| es - Sellador mineral bioactivo de conductos radiculares..... | 6 - 7 |
| de - Bioaktiver mineralischer Wurzelkanalsealer | 8 - 9 |
| it - Sigillante minerale bioattivo per canali radicolari..... | 10 - 11 |
| pt - Selante de canal radicular mineral bioativo | 12 - 13 |
| nl - Bioactieve minerale afdichting van het wortelkanaal | 14 - 15 |
| pl - Biologicznie czynny mineralny uszczelniacz kanałów korzeniowych..... | 16 - 17 |
| el - Βιοδραστικό ανόργανο υλικό έμφραξης ριζικού σωλήνα | 18 - 19 |
| bg - Биоактивен минерален уплътнител за коренови канали | 20 - 21 |
| cs - Bioaktivní minerální těsnicí hmota pro kořenové kanálky | 22 - 23 |
| da - Rodkanalforsegler med bioaktive mineraler..... | 24 - 25 |
| et - Bioaktiivne mineraalne juurekanali täidis | 26 - 27 |
| fi - Bioaktiivinen juurikanavan mineraalitiivistysaine | 28 - 29 |
| hu - Bioaktív ásványi gyökércsatorna tömítőanyag..... | 30 - 31 |
| lt - Biologiškai aktyvus mineralinis šaknies kanalo sandariklis..... | 32 - 33 |
| lv - Bioaktīvs minerālu saknes kanāla blīvētājs..... | 34 - 35 |
| ro - Sigilant bioactiv de canal radicular pe bază de minerale | 36 - 37 |
| sk - Bioaktívna minerálna výplň koreňových kanálikov..... | 38 - 39 |
| sl - Bioaktivno mineralno tesnilo za koreninske kanale | 40 - 41 |
| sv - Bioaktiv mineralbaserad sealer för rotkanal..... | 42 - 43 |
| hr - Bioaktivno mineralno sredstvo za brtvljenje korijenskih kanala | 44 - 45 |
| no - Bioaktiv mineralrotkanalforsegling | 46 - 47 |

BioRoot™ Flow

Active Biosilicate Technology™

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Bioactive mineral root canal sealer

PRESENTATION

| BioRoot™ Flow - REF 340S/2 | BioRoot™ Flow – REF 340S/0.5 |
|--|--|
| 1 syringe of 2.0g of paste contained in a foil | 1 syringe of 0.5g of paste contained in a foil |
| 20 tips contained in a sachet | 5 tips contained in a sachet |
| 1 finger-grip contained in a sachet | 1 finger-grip contained in a sachet |

COMPOSITION

White to slightly grey paste containing tricalcium silicate, propylene glycol, povidone, calcium carbonate, AEROSIL, zirconium oxide, acrylamide / sodium acryloyldimethyltaurate copolymer, isohexadecane, polysorbate.

INTENDED USE

Permanent root canal sealer.

INDICATIONS FOR USE

Permanent root canal filling in combination with gutta-percha points in case of inflamed or necrotic pulp in the course of initial treatment or retreatment procedure

TARGET GROUP(S)

Target population

Mature permanent teeth.

Intended users

Dental healthcare professional.

CLINICAL BENEFIT

Preservation of the tooth functionality

Summary of safety and clinical performance is available, using Basic UDI-DI Code 37601407000005GB, in the European database on medical devices (Eudamed): <https://ec.europa.eu/tools/eudamed> or upon request (sscp@septodont.com) until EUDAMED will be functional.

PERFORMANCES CHARACTERISTICS

Technical performance

Hermetically seal the root canal.

Technical specifications

| | |
|---------------|---------------------------------|
| Working time* | 1 hour ($\varnothing > 17$ mm) |
| Setting time* | 45min – 8h |
| Radiopacity* | ≥ 3mm Al |

*Performed according to the EN ISO 6876:2012

Clinical performance

Permanent root canal obturation of mature permanent teeth.

CONTRA INDICATION

- Hypersensitivity to any component.
- Immature or temporary teeth.

EXPECTED LIFETIME

At least 5 years.

INSTRUCTION FOR USE

Placement into the root canal

1. Perform the root canal cleaning and shaping procedure using standard endodontic procedures.
2. Take out the syringe from the protective foil (set aside the foil for further storage).
3. Remove the syringe cap from the syringe. Place the tip by twisting it into the threads of the Syringe Tip to “lock” it. Tips are flexible and can be bent to facilitate the access to the root canal.
Precaution: Ensure that the tip is correctly locked to avoid tip disassembly.
4. The finger grip can be locked to the syringe to facilitate the extraction of the product.
Precaution: Ensure that the finger-grip is correctly locked to avoid finger-grip disassembly.
5. Check the free flow of the cement by extruding a small amount of cement before the procedure. In case of clogging, the syringe should be discarded.
6. Insert the tip of the syringe into the canal, pull back 2mm and inject the required quantity of BioRoot™ Flow.
7. Complete the obturation by inserting the gutta-percha master cone (single cone technique) or several gutta-percha points (lateral condensation technique).
8. It is recommended to take X-rays to ensure that root canals have been thoroughly filled and that the working length does not exceed 2 mm from the apex.

Cleaning and storage of the syringe

1. Unscrew the tip and dispose of it.
2. Remove any excess of paste from the syringe with a dry and clean tissue and recap the syringe.
Precaution: Ensure that the cap is correctly locked to avoid cap disassembly.
3. Clean the syringe and the finger-grip with ethanol and store them in the protective foil.

Retreatment procedure

Remove the existing sealer (BioRoot™ Flow or other sealer) using conventional removal techniques. Then follow instructions for use given in the section “Placement into the root canal”.

WARNINGS AND PRECAUTIONS FOR USE

Precautions for use

- Only use the supplied tips and finger-grip with BioRoot™ Flow, to avoid component disassembly during the procedure.
- Use the paste within 4 months after the first opening.
- If the syringe, the cap, or the finger-grip falls on the floor clean it with ethanol to avoid contamination.
- BioRoot™ Flow may cause eye and skin irritation. Do not use BioRoot™ Flow without appropriate personal protective equipment.
- Use a rubber dam to ensure that the cement is not in contact with the mucosa or swallowed.
- In order to avoid performance issues (paste extrusion, consistency of the paste) during the next endodontic procedure, recap the syringe and store in the foil after use.

Warnings

- If the tip falls on the floor dispose it to avoid contamination.
- Do not reuse the extruded paste on several patients to avoid cross contamination.
- The working length in the root canal should not exceed 2 mm from the radiological apex to avoid apical extrusion.
- Using BioRoot™ Flow without gutta-percha points jeopardizes endodontic and retreatment procedures.
- The tip is intended for single-use. Discard the disposable tips

- after each application to avoid any risk of cross contamination.
- Clean the syringe and the finger-grip after use to avoid cross contamination.
 - If a post placement is required, it should be performed after complete hardening of BioRoot™ Flow, i.e. after 7 days.

SIDE-EFFECTS

Possible hypersensitivity to any component of BioRoot™ Flow. Please report any serious incident that occurred in relation to the device to Septodont (vigilance@septodont.com) and to the local competent authority for medical devices. Please report all other incidents to Septodont (vigilance@septodont.com).

DISPOSAL

Dispose the used tips in a biohazardous waste bin. No specific recommendations for disposal of other BioRoot™ Flow components.

STORAGE

Store between 5°C and 30°C; protect from moisture.

ACCESS TO THE INSTRUCTIONS FOR USE

Please check regularly to see if the instructions for use have not been updated.

The latest electronic instruction for use as well as historical versions can be downloaded and printed directly from the website www.septodontcorp.com/bioroot-flow-eifu/.

A paper copy can be requested through the website www.septodontcorp.com/bioroot-flow-eifu/.

SYMBOLS

| | |
|---|---------------------------------|
|  | Medical Device |
|  | Unique device identifier |
|  | Period-after-opening |
|  | Authorised Representative in CH |
|  | UK responsible person |
|  | UK marking |

 Callisto Regulatory Consulting Ltd
5 Garden Court
Lockington Hall
Main Street
Lockington
DE74 2RH
UK

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